

Your Professional Recovery System at Home

The Physiolab S1 Cold-Compression Recovery System allows you to continue the same professional care you received in hospital in your own home.

Specifically designed for **knee, hip, shoulder** and **ankle** surgery patients, our advanced cold-compression system combines continuous controlled cooling therapy with targeted intermittent compression to effectively reduce swelling and alleviate pain.

How **Physiolab S1** Works

Using precisely controlled temperature and pressure, with real time feedback, the Physiolab S1 is able to create the optimum healing conditions as described in the published literature. The result is a cold-compression therapy that is safe, effective and comfortable.



Effective pain management



Swelling reduction & inflammation control



Quicker rehabilitation & recovery



	Optimum Therapeutic Zone	Accurately Controlled Temperature	Controlled Pressure	Perceived Comfort
PHYSIOLAB S1	✓	✓	✓	✓
ICE & WATER	✗	✗	✗	✓
GEL PACK	✗	✗	✗	✓

Data from: "A randomised crossover trial of five cryocompression devices' ability to reduce skin temperature of the knee" (Belsey et al, January 2024) and other research by the University of Winchester (data on file).

Physiolab S1

Quick Start Guide

Easy to follow instructions to get you started

Pre-Set Protocols

Enabling optimal healing conditions at home

Faster Recovery

By reducing pain, swelling & inflammation



PROUD TO
SUPPLY TEAM GB



Physiolab Home Recovery Package

- Physiolab S1 device and therapy pack
- Delivered where you need it
- Customer service for support
- Quick start guide and detailed instructions
- Convenient & cost effective ice solution available
- Pre-set programs for safety & efficacy

DISCOUNT CODE

C-ORTHINN

Please refer to our website for pricing. Use the unique code below at checkout to receive your discount.

To place your order please visit physiolab.com



Share your experience



physiolabtechnologies



physiolabtechnologies



physiolab_UK



A BRITISH
COMPANY